



CAMTEC SERIAL NUMBER LIST - BEDS

Key

Example Serial #:

BB10CZ003

BB: Bed
10: Year of Manufacture (2010)
C: Month of Manufacture (March)
Z: Model 3975-S (May be one or two letters)
003: Third bed of this model manufactured that month

Note: the letters 'B' and 'C' were re-used.

B: 4339-0 was manufactured through Dec 2003

B: reused for 3975-0 Started in September 2005

C: 4339-S was manufactured through March 2004

C: reused 3975-S with Hubble Motor (amp plug - used the CMD control box) was manufactured from Jan 2006 - Feb 2006

The date on the serial number will identify the model bed based on the year of manufacture.

The Hubble Motor was only used for a short while. These beds were identified by a special serial number letter to track any possible concerns with the motors.

Product Model Code	Product Model Number	Description
A	3954-S	Expandable Scale Bed 37"/39"/48"/54"
B	3975-0	Expandable Bed Central Locking Sys
C	3975-S	Expandable Bed Central Locking Sys, Hubble Motor
D	4448-0	Bed 48" Spread Base 5 Function
E	4448-S	Bed 48" Spread Base 5 Function with scale
F	4454-0	Bed 54" Spread Base 5 Function
G	4854-0	Bed 54" Spread Base 3 Function
H	5048-0	Bed 48" Spread Base 2 Function
I	5054-0	Bed 54" Spread Base 2 Function
J	4354-0	Expandable Bed 39"/48"/54"
K	4354-S	Expandable Scale Bed 39"/48"/54"
L	6048-0	Bed 48" Spread Base 3 Function
M	Control Box	Control Module "Triacs"


N	3954-0	Expandable Bed 37"/39"/48"/54"
O	4954-0	Bed 54" Spread Base 4 function
P	Control Box	Control Module "D" 14 relay
Q	3954-K	Expandable Bed 37"/39"/48"/54"
R	4754-0	Bed 54" Spread Base 2 Function
S	4754-C	Bed 54" Spread Base with casters
T	4748-C	Bed 48" Spread Base with casters
U	4748-0	Bed 48" Spread Base 2 Function
V	3980-S	Expandable Bed, DC, Battery Back Up, Central Locking w/Scale
W	3980-0	Expandable Bed, DC, Battery Back Up, Central Locking w/o Scale
X	4848-0	Bed 48" Spread Base 3 Function
Y	4948-0	Bed 48 " Spread Base 4 function
Z	3975-S	Expandable, Central Locking, 1000 lb Capacity, with scale
AA	4850-0	Low Bed
AB	4475-0	Expandable, 5 Function, 1000 lb Capacity
AC	4475-S	Expandable, 5 Function, 1000 lb Capacity, with scale
AD	4480-0	Expandable, 3 Function, 750 lb Capacity
AE	4480-S	Expandable, 3 Function, 750 lb Capacity, with scale Discontinued
AF	3990-0	Not Used
AG	3990-S	Not Used
AH	Power Drive System	Camtec Power Drive System
AI	3609-0	NOA Bed
AJ	3700-S	Med Surge
AK	2000-0	Enclosure, Zenith Bed
AL	Patient Lift	Patient Transfer Lift/Trapeze
AM		Camtec Low Bed, 500 lb Capacity
AN		Camtec Low Bed, 1000 lb Capacity
AO		DC Bed
AP	3600-0	Low Bed, 500 lb, No Expansion, 3 Function
AQ	3648-0	Expandable Low Bed, 850 # Capacity, 5 Function
AR	3605-0	Low Bed, 500 lb, No Expansion, 5 Function
AS	3649-0	Low Bed, 850 lb Capacity, Expandable, 5 Function
AT	TD	Transfer Device

RCUHLB850


AU	PERF	Perfusion Bed (Tall)	RCUHLB850
AV	3600-EC	Camtec Model 3600-0 Bed with Enclosure	
AW	3649-S	Expandable Low Bed, 850 # Capacity, 5 Function, with scale	
AX	APD	Lever Type Power Drive	
AY	3648-0	Expandable Low Bed, 850 # Capacity, 3 Function	
AZ	3975-S	UL Listed Bed	RCUHLB850
BA	3605-S	Low Bed, 500 lb, No Expansion, 5 function w/ Integrated Scale	
BB	PBBU	Portable Battery Back Up	
BC	BBU	Battery Back Up Part # 400116-0	
BD	3605-SS	Low Bed, 500 lb, No Expansion, 5 Function with Fall Alarm Scale	
BE	3649-SS	Expandable Low Bed, 850# Capacity, 5 Function with Fall Alarm Scale	RCUHLB850
BF	3975-SS	3975-S with Exit Alarm	
BG			
BH			
BI			
BJ			
BK			
BL			
BM			
BN			
BO			
BP			
BQ			
BR			
BS			
BT			
BU			
BV			
BW			
BX			
BY			
BZ	3975-S	Expandable, Central Locking, 1000 lb Capacity, with Fall Alarm Scale	
NA	IPLUL-5	UL Control Box	
NB	BURAC-3	3 Function Control Box	

NC	BURAC-2	2 Function Control Box
ND		
NE		
NF		
NG		
NH		
NI		
NJ		
NK		
NL		
NM		
NN		
NO		
NP		
NQ		
NR		
NS		
NT	NTAC-3	3 Function Control Box
NU	NTAC-5	5 Function Control Box
NV		
NW		
NX		
NY		
NZ		

Model 3605-O/S/SS							
Spread Base Bariatric Low Bed, 500lb capacity, 5 function, non-expandable - Identificaion in serial number is "AR" for 3605-O, "BA" for 3605-S, and "BD" for 3605-SS							
Description	Camtec Item #	JRC Item #	Detail	Other / or / Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex
Control Box	904009-M	101-200-201-0359					
Knee Gatch Actuator	904003-0 39001521	101-200-201-0072	BLACK				
Hi/Low Foot Actuator	904005-0 39001519	101-200-201-0005	RED				
Fowler (back) Actuator	904004-0 39001520	101-200-201-0027	GREEN				
Hi/Low Head Actuator	904005-1 39001518	101-200-201-0007	BLUE				
Cordset	904179-0		Hospital GRD, 5-15 Down/C13-Left				
Scale Readout	904194-0	101-200-201-0185	(FG1257) -SS				
Jctn Box	904193-0	101-200-201-2370	(FG1257) -SS				
Scale Readout	903773-0	101-200-201-1177	(FG1250) -S				
Jctn Box	904113-0	101-200-201-0329	(FG1250) -S				
Total Lock Caster	903931-0	101-200-202-0294					
Steer Lock Caster	903932-0	101-200-202-0278					
Side Board Head End LH	904325-0 39001486	101-200-202-1243					
Side Board Head End RH	904326-0 39001488	101-200-202-1244					
Side Board Foot End LH	904324-0 39001483	101-200-202-1239					
Side Board Foot End RH	904323-0 39001484	101-200-202-1238					
Mattress	903390-0 PVCO- 3684FLNA25ZP	101-200-205-2502	36"x84"x5"				
Headboard	904312-0 39001477	101-200-205-0735					
Footboard	904328-0 39001475	Pending - Jim M/E.C.	Scale				
	39001476		Deluxe				



Parts List	Manual
	3605-O Manual - Rev I
	3605-S Manual - Rev I
**See Historical Versions	3605-SS Manual - Rev C



FG1250 SCALE INSTRUCTIONS


SYSTEM OPERATION

FG1257 SCALE INSTRUCTIONS


FG1257 BASIC OPERATIONAL FLOW

Model 3649-O/S/SS - RCUHLB850							
Spread Base Bariatric Low Bed, 850lb capacity, 5 function, expandable -							
Identification in serial number is "AS" for 3649-0, "AW" for 3649-S, and "BE" for 3649-SS							
Description	Camtec Item #	JRC Item #	Detail	Other / or / Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex
Control Box	904009-M	101-200-201-0359					
Knee Gatch Actuator	904003-0 39001521	101-200-201-0072	BLACK				
Hi/Low Foot Actuator	904005-0 39001519	101-200-201-0005	RED				
Fowler (back) Actuator	904004-0 39001520	101-200-201-0027	GREEN				
Hi/Low Head Actuator	904005-1 39001518	101-200-201-0007	BLUE				
Cordset	904179-0		Hospital GRD, 5-15 Down/C13-Left				
Scale Readout	904194-0	101-200-201-0185	(FG1257) -SS				
Jctn Box	904193-0	101-200-201-2370	(FG1257) -SS				
Scale Readout	903773-0	101-200-201-1177	(FG1250) -S				
Jctn Box	904113-0	101-200-201-0329	(FG1250) -S				
Total Lock Caster	903931-0	101-200-202-0294					
Steer Lock Caster	903932-0	101-200-202-0278					
Asm, Adj Side Rail Release	300587-0		bed width exp release				
Side Board Head End LH	904325-0 39001479	101-200-202-1243					
Side Board Head End RH	904326-0 39001478	101-200-202-1244					
Side Board Foot End LH	904324-0 39001481	101-200-202-1239					
Side Board Foot End RH	904323-0 39001480	101-200-202-1238					
Mattress	Pending - Jim M/E.C	Pending - Jim M/E.C					
Headboard	904312-0 39001477	101-200-205-0735					
Footboard	900780-0 39001474	101-200-201-0214	no staff ctrl pending Jim	or	300780-S	39001475	Kit, CAMTEC, Footboard, Scale
	904327-0 39001474	101-200-205-0611	staff ctrl pending Jim	or	300780-S	39001476	Kit, CAMTEC, Footboard, Deluxe
Staff Control	903947-0	101-200-201-1101					
Staff Control Cable	903949-0	101-200-201-0601					
Cable Harness (to load cells)	904112-0	101-200-201-0192					

RCUHB850BED Assembly



Extended Warranty Program - RCUHLB850-s Bed Frame		
Item Description	Standard Price	Extended Price
Preventative Maintenance - Annual inspection of bed frames	\$425.00	No Charge
Per facility up to 5 bed frames included per year		
(Charged in increments of 5)		
Discounted on site call out fee (includes 1 hour of labor); per occurrence	\$375.00	\$281.25
Discounted on site labor (after 1st hour); per hour	\$96.00	\$67.20
15% discount on non-warranty parts		
Additional Year warranty coverage:	\$500.00	No Charge
Steel structural parts covered for 6 years		
Electronics covered for 2 years		




Parts List	Manual
	RCUHLB850 3649-O Manual - Rev A
	RCUHLB850 3649-S with Scale Manual - Rev A
	RCUHLB850 3649-SS with Scale Alarm Manual - Rev A
Camtec Off-Set Trapeze Manual - Rev A	Camtec RCUHLB850 Manual - Rev B


Model 3954-O/S							
Fixed frame, 1000lb capacity, 5 function, expandable - Identification in serial number is "N" for 3954-0 and "A" for 3954-S							
Description	Camtec Item #	JRC Item #	Detail	Other /or/ Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex
Control Box	904009-M	101-200-201-0359			300466-S		CMD, Discontinued
Knee Gatch Actuator	904305-0 39001577	101-200-201-0221	BLACK	and	904365-0	101-200-201-2371	BLACK F-Molex to M AMP
Hi/Low Foot Actuator	904304-0 39001576	101-200-201-0164	RED	and	904364-0	101-200-201-0137	RED
Fowler (Back) Actuator	904303-0 39001572	101-200-201-0162	GREEN	and	904362-0	101-200-201-0570	GREEN
Hi/Low Head Actuator	904302-0 39001575	101-200-201-0163	BLUE	and	904363-0	101-200-201-0571	BLUE
Scale Readout	901232-0	101-200-201-1166	No longer Available				
Total Lock Caster	903937-0	101-200-202-0297					
Steer Lock Caster	903938-M	101-200-202-0279					
Asm, Adj Side Rail Release	300587-0		bed width exp release				
Side Rail	300309-0	101-200-202-1269	wing type				
	404003-4 39001482	101-200-205-1261	Pad, Navy				
	301012-0		Replacement/Kit, white				
Mattress		550337-CCBA	Bari w/ Fire Barrier 54"x80"x7"		490040-0	550337-CC	37"/48"x84"x7"
		B1000-4884NZBA	Barimatt 1000 Plus 48"x84"x7"				37"/54"x84"x7"
			Comfort Care Mattress:		490036-0	B1000-4884NZBA	36/53"x84"x7" (1k lbs)
Headboard	300242-0 39001477	101-200-205-0733					
Footboard	106689-0	Pending - Jim M.					

Parts List	Manual
	 3954-S expandable 1000lb bed.pdf

Model 3975-0/S/SS							
Fixed frame, 1000lb capacity, 5 function, expands to 48" & 54"x84"L - Identification in serial number is "B" for 3975-0 and "Z" for 3975-S							
Description	Camtec Item #	JRC Item #	Detail	Other / or / Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex
Control Box	904009-M	101-200-201-0359	NU (NAX9265-3)				
Knee Gatch Actuator	904365-0 39001577	101-200-201-0221	BLACK	and	904365-0	101-200-201-2371	BLACK F-Molex to M AMP
Hi/Low Foot Actuator	904364-0 39001576	101-200-201-0164	RED	and	904364-0	101-200-201-0137	RED
Fowler (Back) Actuator	904363-0 39001572	101-200-201-0162	GREEN	and	904362-0	101-200-201-0570	GREEN
Hi/Low Head Actuator	904362-0 39001575	101-200-201-0163	BLUE	and	904363-0	101-200-201-0571	BLUE
Cordset	904179-0		Hospital GRD, 5-15 Down/C13-Left				
Scale Readout	903773-0	101-200-201-1177	FG1250	or	901232-0	101-200-201-1166	used before 08/10/06
Jctn Box	903787-0	101-200-201-0350	for FG1250 Scale				
Total Lock Caster	903633-0 39001560	101-200-202-0296	5"				
Steer Lock Caster	903634-0 39001561	101-200-202-0280	5"				
Asm, Adj Side Rail Release	300587-0		bed width exp release				
Side Rail	300309-0	101-200-202-1269	wing type				
	404003-4 39001482	101-200-205-1261	Pad, Navy				
	301012-0		Replacement/Kit, white				
Mattress		550337-CCBA	Bari w/ Fire Barrier 54"x80"x7"		490040-0	550337-CC	37"/48"x84"x7"
		B1000-4884NZBA	Barimatt 1000 Plus 48"x84"x7"				
			Comfort Care Mattress:		490036-0	B1000-4884NZBA	36/53"x84"x7" (1k lbs)
Headboard	300242-0 39001477	101-200-205-0733	Plastic	or	300365-0	101-200-205-0740	Oak Veneer
	39001514		39" with Tubes				Melamine
	39001516		39" with Bumpers				Melamine
Footboard	39001491	101-200-205-1416	39" Plastic/no scale	or	300366-0	101-200-205-0599	Oak Veneer
	300328-0 39001492	101-200-205-0600	39" For rotating scale				
	39001489		39"				Melamine
Touch-up Paint	908000-0		Spray Paint, Gray 12oz				



Parts List	Manual
	3975-0 Manual - Rev K
	3975-S Manual - Rev K
	3975-SS Manual - Rev E



Model 4339-0/S

Fixed frame, 1000lb capacity, 5 function, non-expandable -
Identification in serial number is "B" for 4339-0 and "C" for 4339-S (2005 and before)

Description	Camtec Item #	JRC Item #	Detail	Other /or/ Need	Camtec Item #	JRC Item #	Detail	
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex	
Control Box	904009-M	101-200-201-0359	(NU) NAX9265-3	or	300466-S		CMD, Discontinued	
Knee Gatch Actuator	904305-0 39001577	101-200-201-0221	BLACK	and	904365-0	101-200-201-2371	BLACK F Molex to M Amp	A D A P T E R S
Hi/Low Foot Actuator	904304-0 39001576	101-200-201-0164	RED	and	904364-0	101-200-201-0137	RED	
Fowler (Back) Actuator	904303-0 39001572	101-200-201-0162	GREEN	and	904362-0	101-200-201-0570	GREEN	
Hi/Low Head Actuator	904302-0 39001575	101-200-201-0163	BLUE	and	904363-0	101-200-201-0571	BLUE	
Scale Readout	903773-0	101-200-201-1166	FG1250	or	901232-0	101-200-201-1166	used before 08/10/06	
Jctn Box	903787-0	101-200-201-0350	for FG1250 Scale					
Total Lock Caster	902021-0	101-200-202-0121	5"					
Steer Lock Caster	902022-M	101-200-202-0001	5"					
Side Rail	300309-0	101-200-202-1269	wing type					
	404003-4 39001482	101-200-205-1261	Pad, Navy					
	301012-0		Replacement/Kit, white					
Mattress		550337-CCBA	Bari w/ Fire Barrier 54"x80"x7"		490040-0	550337-CC	37"/48"x84"x7"	
		B1000-4884NZBA	Barimatt 1000 Plus 48"x84"x7"					
			Comfort Care Mattress:		490036-0	B1000-4884NZBA	36/53"x84"x7" (1k lbs)	
Headboard	903160-1	101-200-205-2494	Oak veneer				903288-0	
Footboard	903161-1 39001490	101-200-205-2495	Oak veneer					

Model 4354-O/S

Fixed frame, 1000lb capacity, 5 function, expandable -
Identification in serial number is "I" for 4354-O and "K" for 4354-S

Description	Camtec Item #	JRC Item #	Detail	Other / or / Need	Camtec Item #	JRC Item #	Detail	
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex	
Control Box	904009-M	101-200-201-0359	(NU) NAX9265-3	or	300466-S		CMD, Discontinued	
Knee Gatch Actuator	904305-0 39001525	101-200-201-0221	BLACK	and	904365-0	101-200-201-2371	BLACK	A D A P T E R S
Hi/Low Foot Actuator	904304-0 39001524	101-200-201-0164	RED	and	904364-0	101-200-201-0137	RED	
Fowler (Back) Actuator	904303-0 39001523	101-200-201-0162	GREEN	and	904362-0	101-200-201-0570	GREEN	
Hi/Low Head Actuator	904302-0 39001522	101-200-201-0165	BLUE	and	904363-0	101-200-201-0571	BLUE	
Scale Readout	901232-0	101-200-201-1166						
Total Lock Caster	903937-0	101-200-202-0297	5"					
Steer Lock Caster	903938-0	101-200-202-0279	5"					
Asm, Adj Side Rail Release	300587-0		bed width exp release					
Side Rail	300309-0	101-200-202-1269	wing type					
Side Rail Pad	404003-4	101-200-205-1261	navy					
Mattress	490034-0	Pending - Jim M.	39" x 48"	or	490035-0	Pending - Jim M.	39" x 54"	
Headboard	903160-1	101-200-205-2494					903288-0	
Footboard	903161-1	101-200-205-2495						


Model 4448-O/S

it frame, 1000lb capacity, 5 function, non-expandable -
Identification in serial number is "D" for 4448-0 and "E" for 4448-S

Description	Camtec Item #	JRC Item #	Detail	Other /or/ Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex
Control Box	904009-M	101-200-201-0359	(NU) NAX9265-3	or	300466-S		CMD, Discontinued
Knee Gatch Actuator	903802-0 39001525	101-200-201-0071	BLACK	and	904365-0	101-200-201-2371	BLACK
Hi/Low Foot Actuator	903803-0 39001524	101-200-201-0061	RED	and	904364-0	101-200-201-0137	RED
Fowler (Back) Actuator	903841-0 39001523	101-200-201-0022	GREEN	and	904362-0	101-200-201-0570	GREEN
Hi/Low Head Actuator	903842-0 39001522	101-200-201-0319	BLUE	and	904363-0	101-200-201-0571	BLUE
Scale Readout	901232-0	101-200-201-1166					
Total Lock Caster	902021-0	101-200-202-0121	5"				
Steer Lock Caster	902022-M	101-200-202-0001	5"				
Side Rail	305007-0	Pending - Jim M/E.C.	Slide type				
Headboard	903162-1	101-200-205-0738			104605-0	101-200-205-0062	Headboard Post
Footboard	903163-1	101-200-205-0594			104607-0	Pending Jim M/E.C.	Footboard Post

A
D
A
P
T
E
R
S

Model 4475-0/S (RC1000)							
Split Frame, 1000lb capacity, 5 function, expands to 48" & 54" x 84"L - Identification in serial number is "AB" for 4475-0 and "AC" for 4475-S							
Description	Camtec Item #	JRC Item #	Detail	Other / or / Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female	or	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex
Control Box	904009-M (NU)	101-200-201-0359	NU	or	300466-S 300466-N*	101-200-201-0111	(M) or N(NU), Used before 2/2010 *Kit with adapters
Knee Gatch Actuator	904305-0 39001577	101-200-201-0071	BLACK	and	904365-0	101-200-201-2371	BLACK
Hi/Low Foot Actuator	904304-0 39001576	101-200-201-0061	RED	and	904364-0	101-200-201-0137	RED
Fowler (Back) Actuator	904302-0 39001575	101-200-201-0022	GREEN	and	904362-0	101-200-201-0570	GREEN
Hi/Low Head Actuator	904303-0 39001578	101-200-201-0319	BLUE	and	904363-0	101-200-201-0571	BLUE
Cordset	904179-0		Hospital GRD, 5-15 Down/C13-Left				
Scale Readout (FG1250)	903773-0	101-200-201-1177		or	901232-0*	101-200-201-1166	*used before 8/10/2006, obsolete
Jctn Box for FG1250 Scale	903787-0	101-200-201-0350					
Total Lock Caster	903937-0*	101-200-202-0297	5" Total Lock				*used before 2/2008
Steer Lock Caster	903938-0*	101-200-202-0279	5" Steer Lock				*used before 2/2008
Asm, Adj Side Rail Release	300587-0		bed width exp release				
Side Rail	300309-0	101-200-202-1269	wing type				
	404093-4 39001482	101-200-205-1261	Pad, Navy	NEW			
	301012-0		Replacement / Kit, white				
Mattress		550337-CCBA	Bari w/ Fire Barrier	550337-CC			54"x80"x7"
		B1000-4884NZBA	Barimatt 1000 Plus				48"x84"x7"
Headboard	300242-0 39001477	101-200-205-0733	Plastic	or ef	490036-0 300366-0	101-200-205-0740	36/53"x84"x7" (1k lbs) Oak Veneer
	39001514		39" with tubes				Melamine
	39001516		39" with Bumpers				Melamine
	300243-4 39001491	101-200-205-1416	39" Plastic/No Scale	or	300366-0	101-200-205-0599	Oak Veneer
Footboard	300328-0 39001492	101-200-205-0600	39" rotating scale				
	39001489		39"				Melamine



Parts List

[**See Historical](#)

Manual

4475-0 MANUAL, EXPANDABLE BARIATRIC BED

4475-S, MANUAL, EXPANDABLE BARIATRIC BED




Model 4480-0/S							
Expandable Split Base Bariatric Bed, 750lb capacity, 3 function, expandable - Identification in serial number is "AD" for 4480-0 and "AE" for 4480-S							
If customer has a Burac 3 Control Box... The actuators had a different plug style/connector (than the AMP or Molex). Everything else was the same. The would have to replace their control box with another BURAC or switch control box and all actuators. The actuators with this style plug are no longer available.							
Description	Camtec Item #	JRC Item #	Detail	Other /or/ Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	903807-0	101-200-201-1073	upgraded pendant used with NTAC 3 or Burac 3 control boxes	OR	902535-0	101-200-201-1072	CMD Amp plug 3 function Pendant,used with CMD control box and can be used with NTAC 3 or NTAC 5 control boxes via adapter(s)
Hand Pendant (for 4480-S)	903950-0 or 11015942	101-200-201-1058	15-pin Serial Female upgraded 5 function Molex pendant for use with NTAC 5	OR	903284-0*	101-200-201-1056	*used before 2/2010 12-pin Molex 5 function Amp plug
Control Box	904007-0	101-200-201-0355	NTAC 3 control box, for 4480-0	OR	903806-0	101-200-201-0328	for early 4480-0 beds with Burac 3 control box
	300466-N	101-200-201-0111	NTAC 5 control box with adapters, for 4480-S				
Knee Gatch Actuator	903802-0 39001525	101-200-201-0071	BLACK	and	904365-0	101-200-201-2371	Adaptor, F Molex to M Amp, KG, Black
Hi/Low Foot Actuator	903803-0 39001524	101-200-201-0061	RED	and	904364-0	101-200-204-0137	Adaptor, Actuator (Hi/Lo) Ft, CMTC, Red
Fowler (back) Actuator	903811-0 39001523	101-200-201-0022	GREEN	and	904362-0	101-200-201-0570	Adaptor, Actuator, CMTC, Green
Hi/Low Head Actuator	903812-0 39001522	101-200-201-0319	BLUE	and	904363-0	101-200-201-0571	Adaptor, Actuator (Hi/Lo) Ft, CMTC, Blue
Cordset	904179-0		Hospital GRD, 5-15 Down/C13-Left				
Scale Readout	903773-0	101-200-201-1177	for 4480-S				
Jctn Bx for FG Scale	903787-0	101-200-201-0350	for 4480-S				
Total Lock Caster	902021-0	101-200-202-0121					
Steer Lock Caster	902022-M	101-200-202-0001					
Asm, Adj Side Rail Release	300587-0		bed width exp release				
Side Rail	300309-0	101-200-202-1269	wing type				
	404003-4 39001482	101-200-205-1261	Pad, Navy				
	301012-0		Replacement / Kit, white				
Headboard	903795-0 39001511	101-200-205-0764	Kit, CMTC, HeadBoard Melamine, 39"				
Footboard	903796-0 39001489	101-200-205-0606	Kit, CMTC, FootBoard Melamine, 39"				



6 Button 3 Function Pendant

FG1250 SCALE INSTRUCTIONS

SYSTEM OPERATION



Parts List	Manual
**See Historical Versions	4480-0 MANUAL - Rev J

Model 4748-C/O

Split frame, 600lb capacity, 2 function, non-expandable -
Identification in serial number is "T" or "U"

Description	Camtec Item #	JRC Item #	Detail	Other / or / Need	Camtec Item #	JRC Item #	Detail
Hand Pendant	902530-0	101-200-201-1065	Being phased out Verify Stock				
Control Box	902528-0	Pending - Jim M/E.C.					
Knee Gatch Actuator	903802-0 39001525	101-200-201-0071		and	904365-0	101-200-201-2371	BLACK F Molex to M Amp
Fowler (back) Actuator	903811-0 39001523	101-200-201-0022		and	904362-0	101-200-201-0570	GREEN
Total Lock Caster	903639-0	Pending - Jim M/E.C.	x4				
Steer Lock Caster	903640-0	Pending - Jim M/E.C.	x2				
Side Rail	305007-0	Pending - Jim M/E.C.	Slide type				
Headboard	903162-1	101-200-205-0738					
Footboard	903163-1	101-200-205-0594					
Mattress	903260-0	Pending - Jim M/E.C.					

ADAPTERS